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Core Access Surgical Technologies

LevaLap 1.0

LAPAROSCOPIC ACCESS DEVICE

EU REP

MedNet EC-REP GmbH
Borkstrasse 10
48163 Münster
Germany

Core Access Surgical Technologies, Inc.
3495 Piedmont Road
Building 11, Suite 905
Atlanta, GA 30305
USA

www.castsurgical.com

Symbols Glossary

Manufactured in the United States

Manufacturer

European Conformity

Do not resterilize

Do not reuse

Do not use if package is damaged

Use-by date

Medical device

Consult instructions for use

Single sterile barrier system with protective packaging outside

Keep away from sunlight

Keep dry

Temperature limit

Sterilized using irradiation

Catalog number

Quantity

Authorized representative in the European Union

Instructions for use

INDICATIONS FOR USE:

The LevaLap Laparoscopic Access Device is intended for use in the peri-umbilical region of the abdominal wall prior to the establishment of a pneumoperitoneum to aid in the insertion of a Veress needle during gynecologic (pelvic) and general surgical (intra-abdominal) laparoscopic procedures.

INTENDED USE:

The LevaLap Laparoscopic Access Device is intended to raise the abdominal wall above critical organs and vessels by means of negative pressure, thus creating an initial additional peritoneal safe space prior to insufflation in preparation for laparoscopic access.

DEVICE DESCRIPTION:

The LevaLap Laparoscopic Access Device is a dome shaped optically clear device which provides a means of access allowing for rapid creation of the pneumoperitoneum during laparoscopic procedures. The LevaLap Laparoscopic Access Device attaches to standard OR suction, and lifts the abdominal wall, creating additional space for Veress needle introduction. A vacuum bypass plug is positioned in the suction port.

CONTRAINDICATIONS:

- The LevaLap Laparoscopic Access Device should not be used under the following circumstances:
 - In patients who are pregnant.
 - In patients who have had surgery over the chosen access site within the previous 10 days.
 - In patients with confirmed symptomatic abdominal hernia.
 - When use of a Veress needle and/or laparoscopic techniques are contraindicated.
 - In pediatric patients.

WARNINGS:

- FOR SINGLE USE ONLY.
- DO NOT REUSE or RESTERILIZE.
- Do not use the device beyond the expiration date printed on the package label.
- Only for use by a qualified specialist.
- Inspect device and packaging for damage or breakage prior to use. If the device or packaging are damaged do not use the device.
- Do not place excessive stress or force on the device. Excessive force or improper use may result in device failure and/or patient injury.
- Do not use the device on morbidly obese patients with a body mass index of more than 40 kg/m² or more than 100 pounds over ideal body weight.
- Always increase the suction vacuum gradually and do not exceed 60 mmHg (80 mbar or 8 kPa). Increasing the suction vacuum too quickly may cause mild skin irritation.
- Do not use the device for more than 3 minutes to prevent skin bruising.
- Shut off suction once abdominal wall lift is considered adequate by the physician.
- Only for use in laparoscopic patients in whom blind, peri-umbilical passage of a 14G Veress needle for establishment of a pneumoperitoneum is appropriate.

CAUTIONS:

- Federal (US) law restricts this device to sale by or on the order of a physician.
- Always prepare site of application using standard techniques including hair removal. Inadequate preparation may result in inability or failure to maintain the abdominal wall in the lifted position.

LevaLap 1.0

PRECAUTIONS:

- The LevaLap Laparoscopic Access Device is intended for use in the peri-umbilical region of the abdominal wall prior to the establishment of a pneumoperitoneum to aid in the insertion of a Veress needle during gynecologic (pelvic) and general surgical (intra-abdominal) laparoscopic procedures.
- For use only with 14G Veress needles.
- The LevaLap Laparoscopic Access Device has not been demonstrated to reduce the incidence of complications associated with or resulting from blind passage of a Veress needle for establishment of a pneumoperitoneum.
- The LevaLap Laparoscopic Access Device should not be placed on the surface of the abdominal wall skin until the preparatory scrub has dried.

STORAGE:

- The product must be stored at room temperature with normal humidity conditions.
- Keep packaging dry and store out of direct sunlight.

DIRECTIONS FOR USE:

- Using standard sterile technique, remove the LevaLap Laparoscopic Access Device from the sterile packaging and place into the sterile field.
- Prepare the patient's skin surface per standard practice and position the LevaLap Laparoscopic Access Device centrally over the area for intended access.
- Connect an external vacuum system onto the LevaLap Laparoscopic Access Device suction port using appropriate tubing. Gently pull on the tubing to check for secure attachment between the suction port and the vacuum.
- Gently engage the negative pressure source and gradually increase the negative pressure only until the abdominal wall rises to make contact with the septum. Shut off suction once abdominal wall lift is considered adequate by the physician. Do not exceed 60 mmHg (80 mbar or 8 kPa).

NOTE: Keep the LevaLap Laparoscopic Access Device as firmly as possible against the patient's skin using the stabilizing wings during this step to ensure adequate suction is achieved.

- Insert a Veress needle through the LevaLap Laparoscopic Access Device septum.
- Determine that the Veress needle has been inserted to the desired depth and proceed with insufflation per standard technique.
- Gradually release the negative pressure until the initial pneumoperitoneum has been adequately insufflated.
- Once appropriate insufflation has been achieved, remove the Veress needle.
- Turn off the vacuum system and remove LevaLap Laparoscopic Access Device from skin.
- Dispose of used device as per hospital policy.

RESIGNATION OF GUARANTEE:

Core Access Surgical Technologies, Inc. accepts no responsibility for claims deriving from the misuse of the product.

A glossary explaining the symbols used on the labels can be found at www.cast surgical.com.

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